



AimTab™ KETONE TABLETS

**CLIA
WAIVED
FOR
URINE**



AimTab™ Ketone Tablets are for the semi-quantitative determination of ketones (acetoacetic acid and acetone) in urine, serum, plasma and whole blood.

**Ketone Testing
Made Easy!**

Same test procedure as Acetest®, a registered trademark used by Bayer HealthCare LLC.

AimTab™ Ketone Tablets, REF #13100

100 TABLETS
PER BOTTLE



RELATED PRODUCTS

CONTROLS

| | SIZE | REF# |
|------------------------------|-------|-------|
| Ketone Serum Controls | 2x2mL | 13112 |

- ⇒ Negative control (1x2mL)
- ⇒ Positive control (1x2mL)
- ⇒ For use with Ketone tablet tests
- ⇒ Store at 2° to 8°C





| | | |
|---------------------------|-------|-------|
| UA Liquid Controls | 2x5mL | 51116 |
|---------------------------|-------|-------|

- ⇒ Level 1 (1x5mL)
- ⇒ Level 2 (1x5mL)
- ⇒ Used for Ketone control
- ⇒ Store at 2° to 8°C



READING RESULTS¹

| POSITIVE RESULTS = PURPLE COLOR | | NEGATIVE RESULTS = NO PURPLE COLOR | |
|---|--|--|---|
|  | The test is POSITIVE IF A PURPLE COLOR DEVELOPS ON THE TABLET | The test is NEGATIVE IF NO PURPLE COLOR DEVELOPS ON THE TABLET |  |
| READ | | | |
| urine at 30 seconds | | serum & plasma at 2 minutes | whole blood at 10 minutes |



Look at these great features:

- ⇒ Quick & easy ketone results for testing urine, serum, and whole blood.
- ⇒ Detects as little as 5 mg/dL of acetoacetic acid in urine.
- ⇒ Specific for acetoacetic acid and acetone and will not react with betahydroxybutyric acid.
- ⇒ CLIA moderately complex when testing serum, plasma or whole blood.
- ⇒ CAP / CLSI formatted procedures available via fax or email - just add your name!
- ⇒ Store at room temperature.

Sanilabo

SANILABO S.L.

C/San Roque 30, bajo. 46910. Benetusser (Valencia)

☎ 902 100 275 / 963 950 809 📠 902 100 274 / 963 740 484

www.sanilabo.com

Página web

www.sanilaboshop.es

Tienda online

Síguenos en:



ER-0263/2015